



Lytec 2009 Specifications

Workstation System Requirements (recommended):

Equivalent to Pentium 4 3.0Ghz processor or higher, 4GB of available hard disk space, 2GB RAM, 32bit color display, DVD drive.

Workstation Operating system (required):

Windows XP Professional SP2 or higher, Windows Vista Business.

Server System Requirements (recommended):

Equivalent to Pentium 4 3.6Ghz processor or higher, 6GB of available hard drive space, 4GB RAM, 32bit color display, DVD drive.

Server Operating System (required):

Client/Server: Windows 2000 Server or Windows 2003 Server (Enterprise or Standard editions) with latest service packs. Windows 2003 Small Business Server is not supported.

Multi-User/ Professional: Windows XP Professional SP2 or higher, Windows Vista Business.

Network (recommended):

Gigabit full-duplex 10/100/1000 Ethernet router/switch using CAT5e (or CAT6) Ethernet cable. Consult with a licensed network vendor for proper wiring for your business computing environment. Wireless networks are not recommended. If a wireless network is the only option we suggest, at the minimum, an 802.11G wireless network.

Internet (Recommended):

High-speed Cable, DSL, Satellite internet service.

"Windows" User Set-up (Required for Clients who do not utilize CHS for Hardware)

Windows users need full read/write access to the following directories in order to open Lytec, print reports, and use Encoder

C:\Program Files\Lytec 2009

C:\Program Files\Medicode

C:\Documents and Settings\All Users\Application Data\Lytec